







Product data sheet REF 54.35.20

11.12.2025

General data			
article number	54.35.20	Semb lung retr 25cm 20mm fig 1	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826072377
Risk category	1r (EU MDR)	UMDNS-Code	13-649
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,147	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
 			
size	25cm 20mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4301 (X5CrNi18-10) 1.4305 (X8CrNiS18-9)		