







Product data sheet REF 54.35.24

11.12.2025

General data			
article number	54.35.24	Semb lung retr 25cm 24mm fig 2	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826072384
Risk category	1r (EU MDR)	UMDNS-Code	13-649
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,120	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
 			
size	25cm 24mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	1.4305 (X8CrNiS18-9)		