








## Product data sheet REF 55.08.30

11.12.2025

General data			
article number	55.08.30	Bailey aortic rong 12cm 3mm str	
Unit	piece	Customs code (HS code)	90189084
GMDN	35285	GTIN	04046826072681
Risk category	1r (EU MDR)	UMDNS-Code	13-229
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,130	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.53 (catalogue cardio )		
Complete picture		Detail images	
  			
size	12cm 3mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4301 (X5CrNi18-10)		
Others	1.4310 (X10CrNi18-8)		