



## Product data sheet REF 55.10.10

11.12.2025

| General data    |             |                                   |                         |
|-----------------|-------------|-----------------------------------|-------------------------|
| article number  | 55.10.10    | DeBakey vasc dilator Ø 1.0mm 19cm |                         |
| Unit            | piece       | Customs code (HS code)            | 90189084                |
| GMDN            | 11268       | GTIN                              | 04046826072933          |
| Risk category   | 1r (EU MDR) | UMDNS-Code                        | 11-268                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI                       | 4046826;B;R004;01;aj;4W |
| Net weight (kg) | 0,018       | Operating instructions            | GATDR00X                |
| Latex free      | yes         | Preparation instructions          | MPK01-1                 |



|               |                               |
|---------------|-------------------------------|
| Medicon media | 450.00.53 (catalogue cardio ) |
|---------------|-------------------------------|

## Complete picture

## Detail images



|      |              |         |  |
|------|--------------|---------|--|
| size | Ø 1.0mm 19cm | Coating |  |
|------|--------------|---------|--|

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4310 (X10CrNi18-8)              |
| Handle    | 1.4305 (X8CrNiS18-9)              |