








## Product data sheet REF 55.10.30

11.12.2025

General data			
article number	55.10.30	DeBakey vasc dilator Ø 3.0mm 19cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	11268	GTIN	04046826072971
Risk category	1r (EU MDR)	UMDNS-Code	11-268
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;01;aj;4W
Net weight (kg)	0,015	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.53 (catalogue cardio )		
Complete picture		Detail images	
  			
size	Ø 3.0mm 19cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4305 (X8CrNiS18-9)		
Handle	1.4305 (X8CrNiS18-9)		