



Product data sheet REF 55.10.70

11.12.2025

| General data | | | |
|-----------------|-------------|-----------------------------------|-------------------------|
| article number | 55.10.70 | DeBakey vasc dilator Ø 7.0mm 19cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 11268 | GTIN | 04046826073039 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 11-268 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;01;aj;4W |
| Net weight (kg) | 0,027 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK03-1 |



| | |
|---------------|-------------------------------|
| Medicon media | 450.00.53 (catalogue cardio) |
|---------------|-------------------------------|

Complete picture

Detail images



| | | | |
|------|--------------|---------|--|
| size | Ø 7.0mm 19cm | Coating | |
|------|--------------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4305 (X8CrNiS18-9) |
| Handle | 1.4305 (X8CrNiS18-9) |