








## Product data sheet REF 55.11.10

11.12.2025

General data			
article number	55.11.10	Cooley retr 21.5cm serr	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826073053
Risk category	1r (EU MDR)	UMDNS-Code	13-379
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,190	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1
 			
Medicon media	450.00.53 (catalogue cardio )		
Complete picture		Detail images	
  			
size	21.5cm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	3.3535 (AIMG3)		