







Product data sheet REF 55.15.02

11.12.2025

| General data | | | |
|--|---|----------------------------|-------------------------|
| article number | 55.15.02 | DeBakey vasc scrs 25° 16cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 44354 | GTIN | 04046826073282 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-482 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;01;ag;2U |
| Net weight (kg) | 0,029 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 25° 16cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4117 (X38CrMoV15) | | |
| Screw | 1.4305 (X8CrNiS18-9) | | |