



Product data sheet REF 55.15.15

11.12.2025


| General data | | | |
|-----------------|-------------|-----------------------------------|-------------------------|
| article number | 55.15.15 | DeBakey HM vascular scrs 45° 23cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 44354 | GTIN | 04046826073336 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-482 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;01;ag;2U |
| Net weight (kg) | 0,052 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |



| | |
|---------------|-------------------------------|
| Medicon media | 450.00.53 (catalogue cardio) |
|---------------|-------------------------------|

Complete picture

Detail images



| | | | |
|------|----------|---------|--|
| size | 45° 23cm | Coating | |
|------|----------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|--------------------|-----------------------------------|
| Main part Screw | 1.4117 (X38CrMoV15) |
| | 1.4305 (X8CrNiS18-9) |