








## Product data sheet REF 55.30.15

11.12.2025

General data			
article number	55.30.15	DeBakey AT Bulldog clamp 13cm 45° cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	62470	GTIN	04046826073749
Risk category	1r (EU MDR)	UMDNS-Code	10-868
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,018	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.53 (catalogue cardio )		
Complete picture		Detail images	
  			
size	13cm 45°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		