







## Product data sheet REF 55.32.23

11.12.2025

| General data  |  |                                      |   |
|---|--|--------------------------------------|---|
| article number  | 55.32.23   | Ti Cooley AT vasc clamp 11.0cm fig 3 |   |
| Unit  | piece  | Customs code (HS code)               | 90189084  |
| GMDN  | 62468  | GTIN                                 | 04046826450458  |
| Risk category   | 1r (EU MDR)  | UMDNS-Code                           | 15-882  |
| Classification  | 6 (EU MDR)   | Base UDI-DI                          | 4046826;B;R002;04;ah;4K   |
| Net weight (kg)   | 0,014  | Operating instructions               | GATDR00X  |
| Latex free  | yes  | Preparation instructions             |   |
|   |  |                                      |   |
| Medicon media   | 451.04.41 (brochure titanium instruments f Cardiac and Vascular Surgery) |                                      |   |
| Complete picture  |  | Detail images                        |   |
|   |  |                                      |   |
| size  | 11.0cm   | Coating                              | Artificially produced oxide layer (interference colors) on the surface of titanium products |
| Storage conditions  |  |                                      |   |
|    |  |                                      |   |
| Materials   |  |                                      |   |
|   | Material number (DIN designation)  |                                      |   |
| Main part   | 3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))                                      |                                      |   |
| Others  | 3.7035 (ASTM_Grade_2)  |                                      |   |