




Product data sheet REF 55.33.54

25.08.2026

General data			
article number	55.33.54	Neonatal vasc clamp 12.5cm 50° cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826073992
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,019	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
  			
size	12.5cm 50°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		