



Product data sheet REF 55.33.63

11.12.2025

General data			
article number	55.33.63	Castaneda AT Neonatal clamp 14.5cm 20mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826074050
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,025	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
  			
size	14.5cm 20mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		