



Product data sheet REF 55.33.65

11.12.2025

General data			
article number	55.33.65	Castaneda AT Neonatal clamp 14.5cm 25mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826074074
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,025	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1



Medicon media	450.00.53 (catalogue cardio)
---------------	-------------------------------

Complete picture

Detail images



size	14.5cm 25mm	Coating	
------	-------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)