








Product data sheet REF 55.39.15

11.12.2025

General data			
article number	55.39.15	Derra Cooley AT vasc clamp 15mm 17cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826074470
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,038	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
  			
size	15mm 17cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		