







Product data sheet REF 55.40.05

11.12.2025

General data			
article number	55.40.05	DeBakey AT vasc clamp str 22cm serr 55mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826074609
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,065	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
 			
size	22cm 55mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		