







Product data sheet REF 55.44.41

11.12.2025

General data			
article number	55.44.41	Ti DeBakey Glover AT vasc clamp cvd 22.0cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826450472
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,039	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	451.04.41 (brochure titanium instruments f Cardiac and Vascular Surgery)		
Complete picture		Detail images	
			
size	22.0cm	Coating	Artificially produced oxide layer (interference colors) on the surface of titanium products
Storage conditions			
			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		
Others	3.7035 (ASTM_Grade_2)		