



Product data sheet REF 55.59.19

11.12.2025

| General data | | | |
|-----------------|-------------|-----------------------------|-------------------------|
| article number | 55.59.19 | Cooley AT vasc clamp 19.5cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 62468 | GTIN | 04046826074876 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-882 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R002;04;ah;4K |
| Net weight (kg) | 0,042 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |



| | |
|---------------|-------------------------------|
| Medicon media | 450.00.53 (catalogue cardio) |
|---------------|-------------------------------|

Complete picture

Detail images



| | | | |
|------|--------|---------|--|
| size | 19.5cm | Coating | |
|------|--------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13) |