








## Product data sheet REF 55.74.32

11.12.2025

General data			
article number	55.74.32	Satinsky DeBakey AT clamp 32mm 23.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826075125
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,083	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.53 (catalogue cardio )		
Complete picture		Detail images	
  			
size	32mm 23.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		