








Product data sheet REF 55.82.50

11.12.2025

General data			
article number	55.82.50	DeBakey AT vasc clamp 90° 19cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826075286
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,064	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
  			
size	90° 19cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		