






## Product data sheet REF 55.85.15

11.12.2025

| General data  |                                   |                                 |                         |
|---|-----------------------------------|---------------------------------|-------------------------|
| article number  | 55.85.15                          | Dardik AT vasc clamp 15.0cm 15° |                         |
| Unit  | piece                             | Customs code (HS code)          | 90189084                |
| GMDN  | 62468                             | GTIN                            | 04046826075422          |
| Risk category   | 1r (EU MDR)                       | UMDNS-Code                      | 15-882                  |
| Classification  | 6 (EU MDR)                        | Base UDI-DI                     | 4046826;B;R002;04;ah;4K |
| Net weight (kg)   | 0,050                             | Operating instructions          | GATDR00X                |
| Latex free  | yes                               | Preparation instructions        |                         |
|   |                                   |                                 |                         |
| Medicon media   | 450.00.53 (catalogue cardio )     |                                 |                         |
| Complete picture  |                                   | Detail images                   |                         |
|    |                                   |                                 |                         |
| size  | 15.0cm 15°                        | Coating                         |                         |
| Storage conditions  |                                   |                                 |                         |
| No special requirements   |                                   |                                 |                         |
| Materials   |                                   |                                 |                         |
|   | Material number (DIN designation) |                                 |                         |
| Main part   | 1.4021 (X20Cr13)                  |                                 |                         |