



Product data sheet REF 55.86.45

11.12.2025


General data			
article number	55.86.45	Dardik AT vasc clamp 16.0cm 45°	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826075491
Risk category	1r (EU MDR)	UMDNS-Code	15-882
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,050	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1



Medicon media	450.00.53 (catalogue cardio)
---------------	-------------------------------

Complete picture

Detail images



size	16.0cm 45°	Coating	
------	------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)