








Product data sheet REF 56.80.04

11.12.2025

General data			
article number	56.80.04	Mehdorn micro vessel clip 4x1.0mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	17675	GTIN	04046826075910
Risk category	1r (EU MDR)	UMDNS-Code	16-028
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;05;ah;4W
Net weight (kg)	0,001	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
  			
size	4x1.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		