







## Product data sheet REF 56.85.13

11.12.2025

General data			
article number	56.85.13	single clamp 14mm A p/2	
Unit	box/2	Customs code (HS code)	90189084
GMDN	17675	GTIN	04046826076115
Risk category	1r (EU MDR)	UMDNS-Code	16-028
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,001	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.53 (catalogue cardio ) 450.00.60 (catalogue plastic surgery )		
Complete picture		Detail images	
 			
size	14mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		