



## Product data sheet REF 56.85.19

11.12.2025

General data			
article number	56.85.19	single clamp 25mm V p/2	
Unit	box/2	Customs code (HS code)	90189084
GMDN	17675	GTIN	04046826076160
Risk category	1r (EU MDR)	UMDNS-Code	16-028
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,001	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1



Medicon media	450.00.53 (catalogue cardio ) 450.00.60 (catalogue plastic surgery )
---------------	---

## Complete picture

## Detail images



size	25mm	Coating	
------	------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	1.4310 (X10CrNi18-8)