







Product data sheet REF 56.85.47

11.12.2025

General data			
article number	56.85.47	double clamp approximator 8mm A	
Unit	piece	Customs code (HS code)	90189084
GMDN	17675	GTIN	04046826076276
Risk category	1r (EU MDR)	UMDNS-Code	17-675
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,001	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.53 (catalogue cardio) 450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
 			
size	8mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		