



Product data sheet REF 56.85.59

11.12.2025


General data			
article number	56.85.59	double clamp approximator 25mm V	
Unit	piece	Customs code (HS code)	90189084
GMDN	17675	GTIN	04046826076368
Risk category	1r (EU MDR)	UMDNS-Code	17-675
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,001	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1



Medicon media	450.00.53 (catalogue cardio) 450.00.60 (catalogue plastic surgery)
---------------	---

Complete picture

Detail images



size	25mm	Coating	
------	------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4310 (X10CrNi18-8)