







Product data sheet REF 56.86.37

11.12.2025

| General data | | | |
|--|-----------------------------------|--------------------------|-------------------------|
| article number | 56.86.37 | single clip 36mm V p/2 | |
| Unit | box/2 | Customs code (HS code) | 90189084 |
| GMDN | 17675 | GTIN | 04046826076573 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 16-028 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R002;04;ah;4K |
| Net weight (kg) | 0,001 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-1 |
|   | | | |
| Medicon media | 450.00.53 (catalogue cardio) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 36mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4310 (X10CrNi18-8) | | |