






Product data sheet REF 57.01.20

11.12.2025

General data			
article number	57.01.20	Dahlgren cranial rong 21cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826076733
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,250	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media			
Complete picture		Detail images	
			
size	21cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Screw	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)		