



Product data sheet REF 57.04.05

11.12.2025

General data			
article number	57.04.05	Kerrison (Spurling) rong 20cm 5mm 90° up	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826077051
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,280	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
			
Medicon media			
Complete picture		Detail images	
			
size	20cm 5mm 90°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4021 (X20Cr13)		
Others	1.4310 (X10CrNi18-8)		