



Product data sheet REF 57.06.95

11.12.2025

General data			
article number	57.06.95	Kerrison (Schlesinger) rong 18cm 5mm 40° up nf	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826077525
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,245	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
  			
size	18cm 5mm 40°	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4021 (X20Cr13)		
Others	1.4310 (X10CrNi18-8)		