








## Product data sheet REF 57.07.13

11.12.2025

General data			
article number	57.07.13	Kerrison (Micro) rong 18cm 3mm 90° dw nf	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826077617
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,196	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media			
Complete picture		Detail images	
  			
size	18cm 3mm 90°	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4021 (X20Cr13)		
Others	1.4310 (X10CrNi18-8)		