







## Product data sheet REF 57.08.83

11.12.2025

General data			
article number	57.08.83	Kerrison (Micro) rong 20cm 3mm 40° up	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826078393
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,200	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media			
Complete picture		Detail images	
 			
size	20cm 3mm 40°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
<b>Main part</b>	1.4021 (X20Cr13)		
<b>Screw</b>	1.4305 (X8CrNiS18-9)		
<b>Handle</b>	1.4021 (X20Cr13)		
<b>Others</b>	1.4310 (X10CrNi18-8)		