





## Product data sheet REF 57.11.47

11.12.2025

| General data  |  |  |   |
|---|--|--|---|
| article number  | 57.11.47                                 | MC Kerrison (Micro) rong 20cm 2mm 90° dw nf we |   |
| Unit  | piece                                    | Customs code (HS code)                         | 90189084  |
| GMDN  | 35541                                    | GTIN   | 04046826079604  |
| Risk category   | 1r (EU MDR)                              | UMDNS-Code                                     | 13-231  |
| Classification  | 6 (EU MDR)                               | Base UDI-DI                                    | 4046826;B;R003;06;ai;5X   |
| Net weight (kg)   | 0,280                                    | Operating instructions                         | GATDR00X  |
| Latex free  | yes                                      | Preparation instructions                       |   |
|   |  |  |   |
| Medicon media   |  |  |   |
| Complete picture  |  | Detail images                                  |   |
|   |  |  |   |
| size  | 20cm 2mm 90°                             | Coating  | Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors |
| Storage conditions  |  |  |   |
| No special requirements   |  |  |   |
| Materials   |  |  |   |
|   | Material number (DIN designation)        |  |   |
| Main part<br>Screw  | 1.4021 (X20Cr13)<br>1.4301 (X5CrNi18-10) |  |   |