








## Product data sheet REF 57.11.64

11.12.2025

| General data   |  |  |   |
|--|--|--|---|
| article number   | 57.11.64                                 | MC Kerrison rong bay 12cm 4mm 40° up nf we |   |
| Unit   | piece                                    | Customs code (HS code)                     | 90189084  |
| GMDN   | 35541                                    | GTIN                                       | 04046826079710  |
| Risk category  | 1r (EU MDR)                              | UMDNS-Code                                 | 13-231  |
| Classification   | 6 (EU MDR)                               | Base UDI-DI                                | 4046826;B;R003;06;ai;5X   |
| Net weight (kg)  | 0,250                                    | Operating instructions                     | GATDR00X  |
| Latex free   | yes                                      | Preparation instructions                   | MPK20-1   |
|    |  |  |   |
| Medicon media  |  |  |   |
| Complete picture   |  | Detail images                              |   |
|    |  |  |   |
| size   | 12cm 4mm 40°                             | Coating                                    | Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors |
| Storage conditions   |  |  |   |
| No special requirements  |  |  |   |
| Materials  |  |  |   |
|  | Material number (DIN designation)        |  |   |
| Main part<br>Screw   | 1.4021 (X20Cr13)<br>1.4301 (X5CrNi18-10) |  |   |