








Product data sheet REF 57.11.90

11.12.2025

General data			
article number	57.11.90	MC Cushing rong 18cm 2x10mm str	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826079772
Risk category	1r (EU MDR)	UMDNS-Code	15-669
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,250	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media			
Complete picture		Detail images	
  			
size	18cm 2x10mm	Coating	Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4301 (X5CrNi18-10)		
Others	1.4310 (X10CrNi18-8)		