







Product data sheet REF 57.11.91

11.12.2025

| General data | | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|-----------------------------------------|-------------------------------------------------------------------------------------|
| article number | 57.11.91 | MC Cushing rong 18cm 2x10mm 30° ang upw | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 35541 | GTIN | 04046826079789 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-669 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;06;ai;5X |
| Net weight (kg) | 0,250 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-1 |
|   | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 18cm 2x10mm 30° | Coating | Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Screw | 1.4301 (X5CrNi18-10) | | |
| Others | 1.4310 (X10CrNi18-8) | | |