








Product data sheet REF 57.12.02

11.12.2025

General data			
article number	57.12.02	trocar sleeve ti 55mm Ø 11mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45182	GTIN	04046826079802
Risk category	2a (EU MDD)	UMDNS-Code	15-260
Classification	7 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,010	Operating instructions	G571201
Latex free	yes	Preparation instructions	MPK01
 			
Medicon media			
Complete picture		Detail images	
  			
size	55mm Ø 11mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		