









## Product data sheet REF 57.12.75

11.12.2025

General data			
article number	57.12.75	dilator set Ø 18mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45182	GTIN	04046826371890
Risk category	1r (EU MDR)	UMDNS-Code	15-215
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;01;aj;4W
Net weight (kg)	0,100	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media			
Complete picture		Detail images	
  			
size	Ø 18mm	Coating	
Storage conditions			
			
Materials			
	Material number (DIN designation)		
Main part	1.4305 (X8CrNiS18-9)		