







## Product data sheet REF 57.13.01

11.12.2025

General data			
article number	57.13.01	trocar sleeve ti 45mm Ø 14mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45182	GTIN	04046826079888
Risk category	2a (EU MDD)	UMDNS-Code	15-260
Classification	7 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,010	Operating instructions	G571201
Latex free	yes	Preparation instructions	MPK01
			
Medicon media			
Complete picture		Detail images	
  			
size	45mm Ø 14mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		