







Product data sheet REF 57.13.03

11.12.2025

| General data | | | |
|--|-------------------------------------|------------------------------|----------------|
| article number | 57.13.03 | trocar sleeve ti 65mm Ø 14mm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45182 | GTIN | 04046826079901 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 15-260 |
| Classification | 7 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,010 | Operating instructions | G571201 |
| Latex free | yes | Preparation instructions | MPK01 |
|  | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 65mm Ø 14mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.7165 (ASTM_Grade_5_TiAl6V4_(ELI)) | | |