







Product data sheet REF 57.20.62

11.12.2025

General data			
article number	57.20.62	Cushing rong 13cm 2x10mm 30° ang dww	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826081812
Risk category	1r (EU MDR)	UMDNS-Code	15-669
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,060	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media			
Complete picture		Detail images	
 			
size	13cm 2x10mm 30°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4301 (X5CrNi18-10)		
Handle	1.4021 (X20Cr13)		
Others	1.4310 (X10CrNi18-8)		