



Product data sheet REF 57.23.34

11.12.2025

General data			
article number	57.23.34	Caspar rong 16.0cm 4x14mm str fenest serr	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826082369
Risk category	1r (EU MDR)	UMDNS-Code	15-669
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,060	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
			
Medicon media			
Complete picture		Detail images	
			
size	16.0cm 4x14mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4301 (X5CrNi18-10)		
Others	1.4310 (X10CrNi18-8)		