






## Product data sheet REF 57.24.43

11.12.2025

| General data  |                                   |  |                         |
|---|-----------------------------------|--|-------------------------|
| article number  | 57.24.43                          | Cloward ivd rong 15cm 6x10mm 30° ang upw |                         |
| Unit  | piece                             | Customs code (HS code)                   | 90189084                |
| GMDN  | 35541                             | GTIN                                     | 04046826082680          |
| Risk category   | 1r (EU MDR)                       | UMDNS-Code                               | 15-669                  |
| Classification  | 6 (EU MDR)                        | Base UDI-DI                              | 4046826;B;R003;06;ai;5X |
| Net weight (kg)   | 0,069                             | Operating instructions                   | GATDR00X                |
| Latex free  | yes                               | Preparation instructions                 | MPK20-1                 |
|   |                                   |  |                         |
| Medicon media   |                                   |  |                         |
| Complete picture  |                                   | Detail images                            |                         |
|    |                                   |  |                         |
| size  | 15cm 6x10mm 30°                   | Coating                                  |                         |
| Storage conditions  |                                   |  |                         |
| No special requirements   |                                   |  |                         |
| Materials   |                                   |  |                         |
|   | Material number (DIN designation) |  |                         |
| <b>Main part</b>  | 1.4021 (X20Cr13)                  |  |                         |
| <b>Screw</b>  | 1.4301 (X5CrNi18-10)              |  |                         |
| <b>Others</b>   | 1.4310 (X10CrNi18-8)              |  |                         |