








Product data sheet REF 57.26.04

11.12.2025

General data			
article number	57.26.04	MC Ferris Smith rong 30cm 90° dw 4mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826082888
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,320	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	451.04.65 (brochure Thora MIC Thorakosk Mini Inv Surgery)		
Complete picture		Detail images	
  			
size	30cm 4mm 90°	Coating	Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Screw	1.4021 (X20Cr13) 1.4301 (X5CrNi18-10)		