








Product data sheet REF 57.26.53

11.12.2025

General data			
article number	57.26.53	MC Caspar rong 28cm 30° up 4mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826082970
Risk category	1r (EU MDR)	UMDNS-Code	15-669
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,080	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	451.04.65 (brochure Thora MIC Thorakosk Mini Inv Surgery)		
Complete picture		Detail images	
  			
size	28cm 4mm 30°	Coating	Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4301 (X5CrNi18-10)		
Others	1.4310 (X10CrNi18-8)		