







Product data sheet REF 57.28.13

11.12.2025

| General data | | | |
|--|--|--|-------------------------|
| article number | 57.28.13 | Kerrison rong DieBlauen 20cm 3mm 40° up we | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 35541 | GTIN | 04046826404215 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-231 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;06;ai;5X |
| Net weight (kg) | 0,000 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-2 |
|   | | | |
| Medicon media | 450.00.37 (catalogue ENT) 451.04.57 (brochure Medicon DieBlauen dism rongeurs w silicon hdl) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 20cm 3mm 40° | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Screw | 1.4301 (X5CrNi18-10) | | |
| Handle | ELASTOSIL LR3003 (Silikonkautschuk) | | |