






## Product data sheet REF 57.28.14A

11.12.2025

General data			
article number	57.28.14A	Kerrison rong DieBlauen 20cm Comfort 4mm 40° up we	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826460426
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-2
 			
Medicon media	451.04.57 (brochure Medicon DieBlauen dism rongeurs w silicon hdl)		
Complete picture		Detail images	
			
size	20cm 4mm 40°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
<b>Main part</b>	1.4021 (X20Cr13)		
<b>Screw</b>	1.4301 (X5CrNi18-10)		
<b>Handle</b>	ELASTOSIL LR3003 (Silikonkautschuk)		