







## Product data sheet REF 57.28.73

11.12.2025

General data			
article number	57.28.73	MC Kerrison rong DieBlauen 18cm 3mm 90° up nf we	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826404574
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-2
 			
Medicon media	450.00.37 (catalogue ENT ) 451.04.57 (brochure Medicon DieBlauen dism rongeurs w silicon hdl)		
Complete picture		Detail images	
 			
size	18cm 3mm 90°	Coating	Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4301 (X5CrNi18-10)		
Handle	ELASTOSIL LR3003 (Silikonkautschuk)		