








## Product data sheet REF 57.29.63A

11.12.2025

General data			
article number	57.29.63A	MC Kerrison DieBlauen Comfort 23cm 3mm 40° up nf we	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826463977
Risk category	1r (EU MDR)	UMDNS-Code	13-231
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-2
 			
Medicon media	451.04.57 (brochure Medicon DieBlauen dism rongeurs w silicon hdl)		
Complete picture		Detail images	
  			
size	23cm 3mm 40°	Coating	Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4301 (X5CrNi18-10)		
Handle	ELASTOSIL LR3003 (Silikonkautschuk)		